



R-C1

TENS + EMS STIMULATOR DEVICE

Pain Management Solution

- ✓ Pain Relief - Drug Free Technology
- ✓ Muscle Strengthening
- ✓ Massage



Indications for use

Pain Mode (TENS)

To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, neck, upper extremities (arm) and lower extremities (leg) due to strain from exercise or normal household work activities. It is also intended for symptomatic relief and management of chronic, intractable pain and relief of pain associated with arthritis.

Strength Mode (EMS)

- Relaxation of muscle spasms
- Prevention or retardation of disuse atrophy
- Increasing local blood circulation
- Muscle re-education
- Immediate post-surgical stimulation of calf muscles to prevent venous thrombosis
- Maintaining or increasing range of motion

Programs

- 18 Preset TENS Programs
- 15 Preset EMS Programs
- 3 Preset Massage Programs
- 60 Intensity Levels

Specifications

Power Sources	d.c. 6.0V
Frequency	1Hz~125Hz
Pulse Width	30us~380us
Waveform	Biphasic Square Wave
Output Current	0~60mA (at 10 ohm load)
Output Intensity Levels	0-60 levels
Operating Conditions	5°C~40°C; 15%RH~93%RH; 700hPa~1060hPa
Storage Conditions	-10°C~55°C; 10%RH~95%RH; 700hPa~1060hPa
Service Life of Electrode Pads	10 times



Neck Pain

Shoulder Pain



Knee Pain

Ankle Pain



Back Pain

Hip Pain

V2U HEALTHCARE pte Ltd.

Address : 22 Sin Ming Lane #07-88 MidView City, Singapore 573969
Email : business@v2uhealth.com | Tel : (+65) 85 700055

