

# Neuro-Muscular Stimulators (NMS)

## NeuroTrac®



- \* **Basic Pelvic Toner to promote Contenance.**
- \* 3 Custom programmes plus 9 Preset Pelvic Floor programmes:
  - Pelvic Floor pain relief;
  - Urge, Stress, Frequency Incontinence
  - Lack of Sensitivity
  - Pelvic Floor Work Out
  - Endurance
- \* Lock mode function to measure the home compliance (mA and time used).



- \* **Sophisticated Pelvic Toner** has all the features of NeuroTrac® Contingence device, plus:
  - Advanced multi-phase programmes where treatment is split into time phases with different stimulation parameters.
  - Stronger symmetrical stimulation, with the reduced Pulse Width: 50 – 330 µs.
- \* 3 Custom programmes plus 11 Preset Pelvic Floor programmes.
- \* Lock mode function to measure the home compliance (mA and time used).



- \* **Basic stimulator** for Aerobic and Anaerobic muscle training.
- \* 3 Customisable programmes plus 15 Preset Stimulation Programmes for:
  - Warm up, Capillary, Toning
  - Endurance,
  - Calming the muscle.
  - Resistance force output treatment for activities which require prolonged and high activity muscle force: Cycling, Rowing, middle distance Running etc.
  - Maximum force output treatment to produce maximum muscle force output and muscle bulk: weight lifting, judo, ball games, sprint running etc.
- \* Lock mode function to measure the home compliance (mA and time used).



- \* **Extra versatile stimulator** has all the features of NeuroTrac® Sports device, plus:
  - Stronger symmetrical stimulation, with the reduced Pulse Width: 50 – 330 µs.
  - 4 output sockets. Output Intensity [mA] controlled with A and B buttons. A-buttons control 2 sockets on the left, B-buttons control 2 sockets on the right.
  - 21 Preset Sports's programmes and 3 Customisable programmes.
  - Advanced multi-phase programmes where treatment is split into time phases with different stimulation parameters.

